



PRECISIONARY

REFERENCES FOR NON-HUMAN PRIMATE

All of these **peer-reviewed publications** have used the Compresstome® for successful tissue slices. If you have any questions, please contact us:

www.precisionary.com

E-mail: info@precisionary.com

Phone: 617-683-0586

References

1. Abdelaal HM, Kim HO, Wagstaff R, Sawahata R, Southern PJ, Skinner PJ. Comparison of Vibratome and Compresstome sectioning of fresh primate lymphoid and genital tissues for in situ MHC-tetramer and immunofluorescence staining. *Biol Proced Online*. 2015 Jan 7;17(1):2. PMID: 25657614; PMCID: PMC4318225.
2. Chen X, Wolfe DA, Bindu DS, Zhang M, Taskin N, Goertsen D, Shay TF, Sullivan EE, Huang SF, Ravindra Kumar S, Arokiaraj CM, Plattner VM, Campos LJ, Mich JK, Monet D, Ngo V, Ding X, Omstead V, Weed N, Bishaw Y, Gore BB, Lein ES, Akrami A, Miller C, Levi BP, Keller A, Ting JT, Fox AS, Eroglu C, Gradinaru V. Functional gene delivery to and across brain vasculature of systemic AAVs with endothelial-specific tropism in rodents and broad tropism in primates. *Nat Commun*. 2023 Jun 8;14(1):3345. PMID: 37291094; PMCID: PMC10250345.
3. Haran KP, Hajduczki A, Pampusch MS, Mwakalundwa G, Vargas-Inchaustegui DA, Rakasz EG, Connick E, Berger EA, Skinner PJ. Simian Immunodeficiency Virus (SIV)-Specific Chimeric Antigen Receptor-T Cells Engineered to Target B Cell Follicles and Suppress SIV Replication. *Front Immunol*. 2018 Mar 20;9:492. PMID: 29616024; PMCID: PMC5869724.
4. Kalmbach BE, Hodge RD, Jorstad NL, Owen S, de Frates R, Yanny AM, Dalley R, Mallory M, Graybuck LT, Radaelli C, Keene CD, Gwinn RP, Silbergeld DL, Cobbs C, Ojemann JG, Ko AL, Patel AP, Ellenbogen RG, Bakken TE, Daigle TL, Dee N, Lee BR, McGraw M, Nicovich PR, Smith K, Sorensen SA, Tasic B, Zeng H, Koch C, Lein ES, Ting JT. Signature morpho-electric, transcriptomic, and dendritic properties of human layer 5 neocortical pyramidal neurons. *Neuron*. 2021 Sep 15;109(18):2914-2927.e5. PMID: 34534454; PMCID: PMC8570452.

5. Kurz A, Xu W, Wiegel P, Leukel C, N Baker S. Non-invasive assessment of superficial and deep layer circuits in human motor cortex. *J Physiol.* 2019 Jun;597(12):2975-2991. Epub 2019 May 22. PMID: 31045242; PMCID: PMC6636705.
6. Li S, Folkvord JM, Rakasz EG, Abdelaal HM, Wagstaff RK, Kovacs KJ, Kim HO, Sawahata R, MaWhinney S, Masopust D, Connick E, Skinner PJ. Simian Immunodeficiency Virus-Producing Cells in Follicles Are Partially Suppressed by CD8⁺ Cells In Vivo. *J Virol.* 2016 Nov 28;90(24):11168-11180. PMID: 27707919; PMCID: PMC5126374.
7. Li S, Folkvord JM, Kovacs KJ, Wagstaff RK, Mwakalundwa G, Rendahl AK, Rakasz EG, Connick E, Skinner PJ. Low levels of SIV-specific CD8⁺ T cells in germinal centers characterizes acute SIV infection. *PLoS Pathog.* 2019 Mar 21;15(3):e1007311. PMID: 30897187; PMCID: PMC6445460.
8. Webb GM, Li S, Mwakalundwa G, Folkvord JM, Greene JM, Reed JS, Stanton JJ, Legasse AW, Hobbs T, Martin LD, Park BS, Whitney JB, Jeng EK, Wong HC, Nixon DF, Jones RB, Connick E, Skinner PJ, Sacha JB. The human IL-15 superagonist ALT-803 directs SIV-specific CD8⁺ T cells into B-cell follicles. *Blood Adv.* 2018 Jan 23;2(2):76-84. PMID: 29365313; PMCID: PMC5787870.
9. Webb GM, Molden J, Busman-Sahay K, Abdulhaqq S, Wu HL, Weber WC, Bateman KB, Reed JS, Northrup M, Maier N, Tanaka S, Gao L, Davey B, Carpenter BL, Axthelm MK, Stanton JJ, Smedley J, Greene JM, Safrit JT, Estes JD, Skinner PJ, Sacha JB. The human IL-15 superagonist N-803 promotes migration of virus-specific CD8⁺ T and NK cells to B cell follicles but does not reverse latency in ART-suppressed, SHIV-infected macaques. *PLoS Pathog.* 2020 Mar 12;16(3):e1008339. PMID: 32163523; PMCID: PMC7093032.
10. Xu F, Shen Y, Ding L, Yang CY, Tan H, Wang H, Zhu Q, Xu R, Wu F, Xiao Y, Xu C, Li Q, Su P, Zhang LI, Dong HW, Desimone R, Xu F, Hu X, Lau PM, Bi GQ. High-throughput mapping of a whole rhesus monkey brain at micrometer resolution. *Nat Biotechnol.* 2021 Dec;39(12):1521-1528. Epub 2021 Jul 26. PMID: 34312500.
11. Xu W, Baker SN. In vitro characterization of intrinsic properties and local synaptic inputs to pyramidal neurons in macaque primary motor cortex. *Eur J Neurosci.* 2018 Aug;48(4):2071-2083. Epub 2018 Aug 4. PMID: 30019413; PMCID: PMC6175011.